

Single use Urodynamic Catheters

MMS 5000 line

Latex Free

MMS 5000 urodynamic catheters enable ease of patient catheterisation, simultaneous bladder filling and pressure measurement. Our design ensures accurate pressure measurement and due to catheter diameter greater convenience and patient comfort.

MMS also offers cystometry and UPP sets comprising bladder pressure and filling catheters or bladder filling and rectal catheters in one pack.

The range includes different French sizes; all catheters are fitted with winged luer connections for secure and easy attachment to the transducer or filling line. Each filling catheter in the range has a diamond-shaped eye to enable easy disengagement of the pressure catheter following insertion.

The materials used are medical grade PVC and/or Polyurethane, which is suitable for patients with latex hypersensitivity. All MMS 5000 Cystometry and UPP catheters are sterile.

Benefits of MMS 5000 catheters

Strong but flexible

PVC tubing provides a unique combination of strength and flexibility. This makes introduction of the catheter without kinking very easy.

No extension lines

MMS 5000 catheters are presented with long lines, so no additional extension lines are needed.

User friendly

The cm marking and color coding enable fast procedure setup and easy positioning of the catheter.

Universal Connectors

All MMS 5000 catheters have standard luer lock connectors and are compatible with most urodynamic systems.


Complete range

The MMS 5000 line provides a wide range of catheters in several French sizes suitable for different investigations.

MMS 5000 line


8 Fr Double lumen Cystometry catheter, box of 10

MMS 5702




6 Fr Double lumen Cystometry catheter, box of 10

MMS 5716




8 Fr Double lumen UPP catheter, box of 10

MMS 5704




9 Fr Triple lumen Cystometry/UPP catheter, box of 10

MMS 5331



7 Fr Triple lumen circumferential UPP catheter, box of 10

MMS 5333



4.5 Fr Rectal balloon catheter, with introducer, box of 25

MMS 54151



Color	Description	Connection
Neutral ○	infusion line	pump/perfusion
Blue ●	vesical line	Pves
Yellow ●	urethral line	Pura
Red ●	rectal line	Pabd

Table: all MMS 5000 catheters use the same color code

Medical Measurement Systems b.v.
P.O. Box 580
7500 AN Enschede
The Netherlands
T : + 31 - 53 - 480 37 00
F : + 31 - 53 - 480 37 01
E : info@mmsinternational.com
I : www.mmsinternational.com

Medical Measurement Systems USA, Inc.
51 Washington Street, Suite 400
Dover, NH 03820, USA
T : 800 - 236-9310
T : 603 - 750-0037
F : 603 - 750-3155
E : info@mmsusa.net
I : www.mmsusa.net

MMS distributors can be found at: www.mmsinternational.com

Your official MMS distributor:

